






Get a Head Start with Importers in FCS Express

Don't waste time recreating your experiment plots and gates from BD FACSDiva™ and BD Accuri C6™





In 3 Steps you can use your experiment layout as a starting point for your analysis in FCS Express™ with the FREE BD Accuri C6 and BD FACSDiva Importers

- 1 Import a FACSDiva .xml or BD Accuri .c6 file into FCS Express. You can also choose which worksheets or tabs to import.
- 2 Set the page size and orientation of your imported layout in FCS Express.
- 3 Complete the wizard to finish importing and start reporting.



Gating hierarchies are maintained with the ability to add new gates

- Zoom and scaling from BD Accuri acquisition display is maintained



- Unlimited and flexible histogram overlays
- Fully customize imported plots in FCS Express

FCS Express Professional licenses include the BD Accuri and BD FACSDiva importers. Try importing your data for FREE with a 30-DAY DEMO at: www.denovosoftware.com